

Omnilux clear-U Recommended Treatment Protocol

The clear-U handheld LED device comes with 8 Tripaks of specially-formulated skincare that works in combination with the light to help clear acne lesions.

Follow these simple instructions for clearer skin:

- Cleanse the face using the Preparation tripak or a suitable cleanser
- Apply the Treatment serum from the Treatment tripak to your entire face
- Use Omnilux clear-U as directed
- Apply the Aftercare tripak to the entire face over the Treatment serum.
- Wait at least 48 hours before your next clear-U light treatment



We recommend using Photo Therapeutics' Clearing tripaks, as the efficacy of treatment has been proven in studies using this proprietary skincare product.

Is clear-U safe and painless?

Yes. Subjects in our clinical trials have reported some skin redness after the first 1-2 treatments, but this is only temporary and will fade. No other adverse side effects were reported.



OMNILUX clear-U™

for the effective treatment of mild to moderate acne

PHOTO THERAPEUTICS
world leader in light therapy

Contact info for USA:

2720 Loker Ave West Ste C, Carlsbad, California 92010
Tel: 760-607-0488 | Toll free: 800-743-8150 | Fax: 760-607-0288
E-mail: handheld@phototherapeutics.com

Contact info for ROW:

Unit 3 Bonehill Mews, Lichfield Street, Fazeley, Staffordshire B78 3QU
Sales: +44 (0) 1827 252 830 | Support: +44 (0) 161 925 5610
Facsimile: +44 (0) 1827 287 391
E-mail: salesinfo@phototherapeutics.com

www.phototherapeutics.com

Developed from FDA cleared technology used by dermatologists worldwide

Omnilux clear-U™ has been developed from an FDA cleared medical device that has been used by dermatologists worldwide for the successful treatment of acne.

What is clear-U?

Omnilux clear-U uses light technology in the form of light emitting diodes to:

- Destroy the bacteria responsible for acne
- Stimulate healthy tissue growth
- Decrease inflammation and pore size
- Limit oil production in the future

When should I use clear-U?

Use Omnilux clear-U at the first sign of an acne breakout. The treatment program consists of 8 treatment sessions spread over 4 weeks, 2 treatments a week. Use the blue light on the first treatment session; wait 2 days, then use the red light. Repeat this for the next 3 weeks.

Who can benefit from clear-U?

Clear-U benefits anyone suffering with mild to moderate inflammatory acne.

How does the treatment work?

One of the key bacteria responsible for acne is Propionibacterium acnes (P.acnes). This bacteria produces natural chemicals called porphyrins, which are sensitive to light at specific wavelengths. When your skin is flooded with blue and red light, the bacteria self destructs and inflammation is decreased. Your acne fades and your skin heals. It's as simple as that!

OMNILUX clear-U™



ACTUAL SIZE:
5"H x 2.5"W x 1.5"D, 2.9oz.



Photos courtesy of Lee Sy, Park KH,
Department of Dermatology, NMC, Korea

"The Clear-U handheld unit provides our patients with the tried and true results of LED technology in the palm of your hand. Finally an at-home unit that really works to treat acne! LED technology is safe, effective and easy to use. What more could you ask for?"

Glynis Ablon, M.D., F.A.A.D., Assistant Professor UCLA,
Director Ablon Skin Institute

What results can be expected from a series of treatments?

Best results will be seen between 4-8 weeks after treatment. The light stimulates natural processes that occur after the treatment has stopped. Studies have shown complete clearance of lesions, however, on average you should expect 7 out of 10 acne lesions to clear.

How many treatments will I need?

The course of treatments is 2x20 minute treatments per week for 4 weeks of alternating blue and red light. However, because the light will not damage your skin, there is no limit on the number of treatments you can have.

How long do the effects last?

For most people, their acne is clear for six months. In some cases, follow-up treatments (typically 2 to 4 treatments) may be necessary.

Is clear-U all I need for my acne?

Omnilux clear-U is not a substitute for good skin hygiene or a doctor's care. We suggest that the light therapy is used in conjunction with a preparation containing salicylic acid. Your skin specialist can discuss this combination treatment with you.

order today...

USA orders: Call 800-743-8150 or email
handheld@phototherapeutics.com

ROW orders: Call +44 (0) 1827 252 830
salesinfo@phototherapeutics.com