

## **See how the Personal Medik8 Titanium Dermastamp can help with your skin problem:**

### **Eye Wrinkles**

The delicate eye area is one of the first places to give away our age and also the most difficult area to treat due to the sensitivity and thinness of the skin.

Piercing the upper layers of the skin with microneedles leads to increased blood flow in the skin and stimulates the synthesis of collagen and elastin. Enhanced penetration of topically applied Medik8 products will also assist in stimulating the production of these skin proteins. Each subsequent treatment gradually reduces the appearance of lines and wrinkles.

Use only a 0.2mm Titanium Dermastamp for areas around the eye due to the thin and delicate nature of the skin in this area.

Notice: To build up skin tolerance to skin needling, allow 3 day intervals between treatments for a minimum of 2 weeks. As your skin becomes adjusted to needling, gradually increase the frequency of use. Do not use more than once a day. If used more than once a day the Medik8 Titanium Dermastamp may cause excessive inflammation and even cause damage.

### **Medik8 Dermastamp Eye wrinkle Protocol**

For the Eye Area

1. **STERILISE:** Disinfect the Titanium Dermastamp by using Medik8 Sanitising Solution directly on the needles before every use and allow to dry.
2. **CLEANSE:** Cleanse the whole eye area by applying Medik8 oneCleanse or Red Alert Cleanse or your preferred Medik8 cleanser, following product instructions. Ensure the skin is left clean and dry.
3. **STAMP:** Hold the Dermastamp (like a pen) close to the needle head, providing greater control and accuracy of application. The free hand can be used to hold the skin taught surrounding the eye to stretch the lines/ wrinkles.

Rest your elbow on a flat surface to avoid movement of the needles when in contact with your skin.

Hold the Dermastamp directly over the targeted wrinkle. Repeatedly press the needles into the skin 20-25 times ensuring the whole length of the wrinkle receives repetitive needling.

**IMPORTANT TIP:** fully remove the Dermastamp from contact with the skin between each application. Do not rotate/ twist the Dermastamp whilst in contact with the skin.

4. **PREVENT:** Apply Medik8 vitamin C serum; C-Tetra or CE-Tetra directly over the stamped area and allow to fully absorb.
5. **CORRECT:** Apply Medik8 Retinol Eye TR directly over the stamped area and allow to fully absorb.

Optional: Apply Pretox Eyelift or Dark Circles to the under eye area for a complete eye rejuvenation and younger looking appearance.\*

6. FINISH: Clean the Dermastamp by running under hot water and sanitise as before.

Allow to dry before replacing cap.

\*Application of these products should proceed the application of vitamin C and Retinol serums in line with the Medik8 protocol.

Treatment duration: 5 – 10 minutes

Frequency of use: 2 – 7 times per week, depending on skin tolerance.

To build up skin tolerance to Skin Needling, allow 3 day intervals between treatments for a minimum of 2 weeks. Gradually increase frequency of use, as your skin becomes adjusted to needling.

Do not use more than once a day.

Optimal treatment time: For best results, use at night to allow for maximum absorption of the active formulation