

New: Skinrölla



SKINRÖLLA HOME USE

Skinrölla Home Use is a specially designed device that dramatically increases production of the skin's natural components such as collagen and elastin, creating younger looking skin.

At home, Skinrölla stimulates the release of signals that rejuvenate the skin surface for a more youthful appearance, reducing pore size and improving skin tone.

Skinrölla Home Use is ideal to recommend to clients who wish to treat photo-ageing and hyperpigmentation, premature ageing, fine lines and wrinkles and clients wishing to further enhance the results of their in-clinic treatments.

Skinrölla Home Use has been developed to create a desirable and easy to use device for clients to use. Initially you purchase a complete kit that comes with a travel case and a unique two-phase cleaning and sanitising system and then simply replace just the head each time a fresh set of micro-needles is required (usually every 3-6 months).

The device is washed and sanitised with a special spray after each use as normal and then an additional deep cleaning system is included used every 5 applications for additional reassurance.

Skinrölla Home Use is available in sizes 0.2mm and 0.3 mm

Why the Skinrölla:

- CE marked medical device
- Ergonomic non-slip grip
- Safety features to prevent accidental needle sticking

-Medical grade titanium alloy is harder, sharper and more durable to maintain the efficacy of the device throughout its recommended shelf life.

-Unique Replaceable Head System allows the client to safely remove and exchange the device head every 3-6 months to ensure continued effectiveness and affordability.

Skinrölla Professional and Home Use devices provide a complete system to provide professional, quality micro-needling treatments.

For optimum results Skinrölla treatments can be supported with the use of Medik8 homecare such as Ultimate Recovery Cream and Hydr8 B5 to promote healing. To maximise the results of collagen induction and cell renewal, Medik8 CE Thione or CE Tetra in the morning and Retinol 3TR or 6TR in the evening can also be used for skin priming and after care.

If you are concerned about hyper-pigmentation or to minimise post-inflammatory hyperpigmentation; Medik8 White Balance Click is the ideal choice to support this treatment.