

## **Is the Quasar MD PLUS for me?**

### **Can I prevent wrinkles?**

The Quasar MD PLUS is the perfect way to prevent future wrinkles and promote new collagen production.

### **Will Quasar MD PLUS reduce the signs of stress?**

Quasar MD PLUS will help you prevent unwanted wrinkles with light therapy. The Quasar MD PLUS is a medical tool specifically designed to reduce the wrinkles resulting from stress and fight the signs of aging by promoting new collagen production. Our thirties and forties are critical years when prevention is encouraged. The increased collagen will prevent deep wrinkles from forming.

### **Will the Quasar MD PLUS save me money versus chemical treatments, cosmetic surgery, and injectables?**

Absolutely! Save hundreds to thousands of dollars per year versus chemical treatments, cosmetic surgery, and injectables. You can use your Quasar MD PLUS in the comfort of your own home with no downtime or risk of scarring.

## **FAQ's - Red light**

### **What is LED Light Therapy?**

Photo rejuvenation or LED Light Therapy is the use of various colors of light emitting diodes (LEDs) that create an anti-aging effect on the skin. Exact wavelengths of light emitted from LED's stimulate and energize cellular activity, which in turn creates collagen and elastin for Anti-Aging benefits. LED Light Therapy has become recognized as one of the safest, quickest, and most affordable ways to achieve younger-looking skin. It's referred to by several different names including: Photo-Rejuvenation, Anti-aging LED Light Therapy, LED Light Therapy or Red/Yellow Light Therapy.

### **Is LED Light Therapy a laser?**

The LED Light Therapy is not a laser treatment which is why it is safe for home use or professional use. Lasers utilize a highly focused and coherent beam of a specific light frequency, and can generate enormous heat and energy at the focal point. Light from LEDs is not focused, is non-coherent, and generates no heat. LED Light Therapy is non-invasive, with no down time, requires no surgery, no harsh chemicals, no lasers, no numbing cream, or recovery time.

LED Light therapy works, it's affordable and you can use your spa-quality unit in the privacy of your own home or in your health care facility. And it takes just a few minutes a day to use. LED Light Therapy is safe for all skin types and ages. It is today's natural skin care at its finest.

### **Is there any scientific proof of Photo Rejuvenation or LED Light Therapy?**

Unlike many technologies on the market today, LED Light Therapy has been thoroughly tested, evaluated and validated by hundreds of universities, medical labs and most importantly NASA. There have actually been over 2000 reports and documents written on this technology.

It has been proven to stimulate the production of collagen and elastin, smooth skin texture, improve skin laxity, increase circulation and lymph drainage and reduce the appearance of fine lines, wrinkles and superficial hyperpigmentation.

- Minimises fine lines and wrinkles
- Skin is smoother, firmer and plump
- Evens out skin tone
- Increases hydration
- Soothes irritation and redness
- Stimulates blood circulation
- Accelerates the healing process
- Restores natural luminosity and glow

### **Is Photo Rejuvenation or LED Light Therapy safe?**

Yes, it is safe and completely non-invasive. It has been used for years in numerous medical procedures, medical spas, the offices of aestheticians and dermatologists, and by consumers around the world.

### **Are there any side effects or downtime involved from the treatments?**

No. There have been no reported adverse side effects like redness or irritation from the treatment, and the light does not damage skin in any way.

### **Will LED treatments hurt or be painful?**

Not at all. Because there's no heat involved, there's no pain! Unlike lasers or other ablative (skin harming) treatments, LED treatments are non-thermal, non-ablative, and non-invasive. Clients report no discomfort or pain associated with LED treatments.

### **Is LED Light Therapy safe to use around my eyes?**

LEDs are very bright and therefore it is NOT advisable to look or stare directly into the light. No studies have been found showing any harmful effects to the eyes when using our product correctly. Keep the eyes closed during use on the face or use the comfort goggles enclosed and the light will not cause damage to the eye.

### **How do I prepare my skin before and after using the light device?**

Before use, simply cleanse your skin as usual. Do not apply products prior to use. Use the unit on clean skin only. After treatment you may apply your favourite antioxidant or peptide-rich moisturiser. It is also beneficial to apply daily sun screen to protect your investment.

### **How long until I see results?**

Everyone's skin varies, but most people begin to feel and see a difference after ten days to 2 weeks of daily use. Often the first improvement noticed is that the skin pores are becoming smaller, and for ladies, their makeup starts going on much smoother. Because of differences in damaged skin, from age, skin condition, diet, alcohol consumption, smoking habits, etc., some people will see the benefits almost immediately and others will need more treatments. For collagen production and anti-aging, you should see results in 3 to 4 weeks. However, some individuals experience positive results in a shorter timeframe. Skin texture and tone improve from the first treatment. Pigmentation usually fades over the course of 4 to 5 weeks or more of treatments.

The most important things to remember are is that the Facial Secret technology works, so it's just a matter of consistent daily use to determine 'How long will it take?' Also it's very important to take a close-up picture of your face before you start using your Facial Secret so that you can see the improvements clearly in the weeks to come by comparing your face in the mirror to the picture.

### **Do LED Treatments work on everyone?**

LED treatments are proven to work on all skin types. There is an over 90% success rate for people using LED light therapy.

### **Can I combine LED treatments with other skin care products?**

Yes, you can. LED treatments can be used alone or with other skin care treatments. In fact your results will happen faster and be more dramatic if used with a high quality anti-aging or moisturizing skin care product. Many physicians have found LED treatments to complement other skin care treatments.