

Quasar

QUICK START GUIDE

1. First cleanse the skin in the area you wish to treat.
2. Plug 12v power supply connector into the opening at the base of the handle. Then plug the power transformer into a wall socket.
3. Push the Red button firmly to turn the Quasar on.
4. Treat your first area for 3 minutes. Place the light gently against the skin and hold for 3 minutes. (The Quasars have a timer. The Quasar will sound a tone at 3 minutes and it will shut off. This is a safety feature.)
5. Push the Red button firmly to turn the Quasar back on and treat your next area repeating step 4 until all areas have been treated.
6. **Recommended treatment is 5 days per week, treating each area for 3 minutes.**

Quasar

NOTE: Your Quasar comes with an optional clear protective lens. YOU MAY leave the Protective lens on your Quasar during treatments to keep your Quasar clean OR YOU MAY use the Quasar without the protective lens.

NOTE: Check that the Quasar is operating correctly by looking at the LEDs from an angle, so that it is not too bright for your eyes. All the LEDs should be lit except 6 infrared LEDs which are not visible to the human eye.

NOTE: Always clean your Quasar after each use. Unplug your Quasar, wipe the tops of the LED bulbs and the surface of the rim with alcohol wipes or a small amount of rubbing alcohol on a cotton pad or cotton ball. Wash the protective lens with warm soapy water. Replace when dry.