




Integrated under bust effect, for the best surgical result and elegantly shaped breasts

Ensures the best surgical result, without compromising user comfort

Ava

SEAMLESS POST-OP BRA



For augmentation

COMPRESSION LEVEL 2 - FIRM



Support/
compression zones



Breathable




4-way stretch




Seamless



Fully openable



High back for optimum support. Support zones at the back and sides.



Fully openable for easy application

The Ava bra supports, sculpts and contours the breasts, thereby ensuring best surgical result and elegantly shaped breasts.

The 4-way stretch material always adapts perfectly to the breast and body contour. This way the Ava bra contributes to a gentler healing process.

With the front opening and openable straps, it is easy to apply the bra for both health care professionals and user, and it allows easy examination of the surgical site.

With the specially designed built-in compression zones (Tytex ZoneTech®), the Ava bra provides an optimum balance between compression and comfort.

Indications:

Medical device for post-surgical use after breast surgery especially breast augmentation. Also suitable for, breast reduction, breast lift and breast reconstruction.

Precautions:

If you have any questions regarding usage and correct fitting, please consult your doctor or nurse. Ensure that any bandages and/or drainage tubes underneath the bra are positioned correctly to avoid pressure points and/or obstructed flow from the drainage tubes. For breast augmentations, a stabilizer band for breast implants may also be required. Consult your doctor or surgeon before use. Always follow the surgeon's recommendations. If any unexpected discomfort or wound irritation occurs, please contact your doctor or nurse.

Recommended sizes

Band Size	Cm	65 cm	70 cm	75 cm	80 cm	85 cm	90 cm	95 cm	100 cm	105 cm	110 cm	115 cm
	In	26 in	28 in	30 in	32 in	34 in	36 in	38 in	40 in	42 in	44 in	46 in
Cup size	A											
	B		S(Cup A-C)	M(Cup A-C)			L(Cup A-C)				XL(Cup A-C)	
	C											
	D											
	E	S(Cup D-F)		M(Cup D-F)			L(Cup D-F)				XL(Cup D-F)	
	F											

* NOTE: Due to possible post surgical swelling and edema, you may need to go up 1-2 sizes compared to your normal bra size. For more compression go down one size and for less compression go up one size. Important: Always measure yourself to find the right size!

Available colours:



Composition:

78 % Polyamide
13 % Elastane
9 % Polyester

Care Handling:



Close zipper before washing. Wash with similar colours in a washing bag. Use of chlorine-based disinfection will decrease durability and function.

Dealer: