

CLEANING AND STORING

1. Wash brush and foam attachments with anti-bacterial soap after use and allow to dry before storing.
2. Do not use alcohol, acetone, petrol, abrasive cleaning agents, etc. to clean the plastic parts.

The foam and brush heads should be replaced regularly, approximately every 3 months (depending on use).

DermaTx Microdermabrasion Creams and replacement attachments are available to purchase separately from authorised stockists or :

WWW.DERMATX.CO.UK



Advanced skincare designed for your skin

“Exfoliation is one of the most important skincare steps for healthy, radiant skin.”

FACIAL TOOL MANUAL



DermaTx[®]
ADVANCED SKINCARE

CLEANSING

CONTAINS:

- 1x Electronic tool
- 1x Cradle
- 1x Exfoliating foam head
- 1x Cleansing brush head
- 1x User manual

BEFORE USE

1. Remove battery cover (pull off bottom section of the tool firmly).
2. Insert 2 AA batteries matching the +/- symbols.
3. Replace the battery cover.

CLEANSING

Cleansing with a powered device has been proven to cleanse up to **4X MORE EFFECTIVELY THAN MANUAL CLEANSING**, helping to gently loosen dirt and oil and clear pores to leave your skin feeling soft and smooth.

Attach **CLEANSING BRUSH HEAD** to the DermaTx Facial Tool. Moisten both the head and the skin with warm water. Apply DermaTx cleanser to your skin or onto the brush head. Switch device on and apply using small circular motions around the face for **1 MINUTE**. Rinse, pat dry and apply a moisturiser.

Use the cleansing brush daily.

DO NOT USE PRIOR TO A DERMATX MICRODERMABRASION TREATMENT.

DERMATX
CLEANSER



CLEANSING
BRUSH HEAD



TIME



www.dermatx.co.uk

MICRODERMABRASION

For use with DermaTx Microdermabrasion Creams.

The massaging action allows the crystals in the cream to gently exfoliate the skin and enables the active ingredients to be absorbed more easily.

MICRODERMABRASION

Cleanse face manually (**NOT WITH A POWERED CLEANSING BRUSH**) before a microdermabrasion treatment and rinse.

Attach the **EXFOLIATING FOAM HEAD** to the device, apply DermaTx Microdermabrasion Cream and dot onto clean damp skin (forehead, cheeks, chin, nose). Switch device on and using circular motions, massage the cream into the skin working different areas of the face for **2 MINUTES IN TOTAL**. Rinse off, dry face and apply a moisturiser.

Perform treatments **1-3 TIMES A WEEK** for optimal results and overall improvement of skin tone and texture.

CAUTION: SEEK ADVICE BEFORE USE IF TAKING ORAL OR TOPICAL RETINOIDS.

DERMATX
CREAM



EXFOLIATING
FOAM HEAD



TIME



www.dermatx.co.uk