

## DERMATX SKIN RESURFACING TOOL

DermaTx Microdermabrasion Creams have been developed for use with the DermaTx Skin Resurfacing Tool which features interchangeable attachments for different applications: microdermabrasion & daily cleansing.

1 x electronic tool - 1 x cradle - 1 x exfoliating foam head - 1 x cleansing brush

### EXFOLIATING FOAM HEAD

For use with DermaTx Microdermabrasion Creams. The massaging action allows the crystals in the creams to gently exfoliate the skin and enables the active ingredients to be absorbed more easily.

### CLEANSING BRUSH HEAD

Cleansing with a powered device has been proven to cleanse up to 4x more effectively than manual cleansing, helping to gently loosen dirt and oil and clear pores to leave your skin feeling soft and smooth.

### BEFORE USE

- 1 - Remove battery cover (pull off bottom section of the tool firmly).
- 2 - Insert 2 AA batteries matching the +/- symbols.
- 3 - Replace the battery cover.

### CLEANING AND STORING

- 1 - Wash brush and foam attachments with anti-bacterial soap after each use and allow to dry before storing.
- 2 - Do not use alcohol, acetone, petrol, abrasive cleaning agents, etc. to clean the plastic parts.

The foam and brush heads should be replaced regularly, approximately every 3 months (depending on use).

DermaTx Microdermabrasion Creams and replacement attachments are available to purchase separately from authorised stockists or:

[www.dermatx.co.uk](http://www.dermatx.co.uk)

Exfoliation is one of the most important skincare steps for healthy, radiant skin

## MICRODERMABRASION



DERMATX CREAM



EXFOLIATING FOAM HEAD



TIME

Cleanse face manually (NOT with the Cleansing Brush) before a microdermabrasion treatment and rinse. Attach the EXFOLIATING FOAM HEAD to the device, apply DermaTx Microdermabrasion Cream and dot onto clean skin (forehead, cheeks, chin, nose). Switch device on and using circular motions, massage the cream into the skin working different areas of the face for 2 minutes in total. Rinse off, dry face and apply a moisturiser. Perform treatments 1-2 times a week (or up to 3x maximum if required) for optimal results.

Caution: seek advice before use if taking oral or topical retinoids.

## CLEANSING



CLEANSER



CLEANSING BRUSH HEAD



TIME

Attach CLEANSING BRUSH HEAD to the device and apply your cleanser. Switch device on and gently brush your skin using small circular motions around the face for 1 minute. Rinse, pat dry and apply a moisturiser.

Use the cleansing brush daily.

Do not use prior to a DermaTx microdermabrasion treatment.

[www.dermatx.co.uk](http://www.dermatx.co.uk)



At-home microdermabrasion treatments - designed for *your* skin

Formulated using the most advanced ingredients.

Immediate results - Achieving smoother, younger looking skin has never been easier!



**DermaTx**  
ADVANCED SKINCARE

