

Precautions:

- Following treatment, the skin will appear flushed. This should subside within a few hours and can be easily covered with makeup.
- Alternatively, only use the Medik8 Titanium Dermastamp at night. - Cease use if you develop a skin rash or irritation and resume treatment after skin feels comfortable again.
- Do not share your Medik8 Titanium Dermastamp with anyone else.
- Do not use broken skin or open wounds.
- Do not use if pregnant.
- Do not use if you are taking blood thinning medication or are a haemophiliac.
- Store out of reach of children.
- Not to be used by anyone under the age of 18.

Contraindications:

Do not use on the following skin conditions:

- Active acne outbreaks or active skin infections.
- Eczema, psoriasis and other chronic conditions.
- Rosacea and blood clotting problems.
- Actinic (solar) keratosis.
- Keloid or raised scars.
- Raised moles or warts.
- Excessively sensitive skin