

**GUARANTEE: 100% Contaminant-Free**

MTS Rollers™ are sterilized with gamma radiation and sealed in a hermetic pouch and are guaranteed to be 100% contaminate-free upon delivery.

#### Cleaning MTS Roller™

Open the cleaning tray and place the white cleaning pad in the larger half of the tray. Fill the cleaning tray with sterilization solution and roll the roller into the pad with medium pressure. Leave the MTS Roller in the cleaning solution for about 15 to 30 minutes – ensure that the roller drum is constantly wet with the cleaning solution. After cleaning, rinse with tap water and return the roller to the protective case and allow it to air dry. Replace the case lid until the next use.



#### Storing MTS Roller™

The MTS Roller is delivered in a protective case and sealed in a sterilized pouch. Tear open the pouch at the notch and discard. After each use (and cleaning) return the MTS Roller™ back to the protective case.

**NOTE:** Always treat the MTS Roller™ gently and avoid any contact with hard objects or surfaces. Replace all rollers that have been dropped. Keep rollers out of reach of children.

#### Important Reminders

Do not use the MTS DermaRoller™ on irritated skin, infected skin, fungal skin infections, active acne, rosacea, eczema, psoriasis, severe solar keratosis, raised moles and/or warts, and/or any open wounds or sores. Patients who are not sure about their skin condition should consult a dermatologist before using the MTS Roller™.