


YOU ARE BEAUTIFUL
- BE ALWAYS



lipoelastic[®]



THE UNDERMENTIONED CERTIFICATES PROVES THE QUALITY OF OUR GARMENTS

ÖKO-TEX

It guarantees the absence of harmful substances in our garments that are yearly tested by the Trial Institute Hohenstein in Germany. The certificate is the guarantee of the quality and the safety for our patients.



ISO 9001- 2001

This international standard is the proof of the effectivity of the management system and the production with the aim to intensify the customers' satisfaction and to fulfill their requests. The observation of the system is yearly supervised by the competent authority.



CE - DECLARATION OF CONFORMITY

LIPOELASTIC® products meet the requirements of the Medical Device Directive 93/42/EEC and are authorized to use the CE mark.

THE UNIQUE COLLECTION OF LIPOELASTIC®

Only LIPOELASTIC® offers the patients the choice of material lines Classic and Micro. Our materials guarantee the ideal compression with the connection to comfort. They guarantee the elasticity and the tensibility to all directions which assures the stable compression and the equal pressure for the operated sections.


CLASSIC



MICRO



MATERIAL COOL EFFECT



The lines are different in material composition. The materials guarantee COOL EFFECT - the optimal way of transporting moisture away from the body.

NEWS BANDAGES DESIGN



The front double underlay is now available in the modern design.



The new trendy lace has a silicone wavelet on the reverse side. Therefore the garment do not throttle and the fixation in this part is still guaranteed.

DESIGN OPTIONS



WITH ZIPPERS

Design with zippers at the sides.



WITHOUT ZIPPERS

Design without any zippers, ideal for the second period of the long-time compressive therapy.



COMFORT

Zippers with hooks underlay for much easier putting on.



VARIANT

Closing by three-position hooks only, which offers the variability of setting the garment after the swelling falling.

The effort of the company LIPOELASTIC® is to satisfy all patients. That's why we produce our garments in 3 colours and sizes of SS through XXXXL. The bandages are produced by the highest quality sewing machines. The sewing technology applied ensures flat, strong and resistant seams. Whole the production includes up to 30 production steps.

HYGIENIC CROTCH TYPES



All male and female garments with trouser-leg



All male and female garments without trouser-leg

CARE:

LIPOELASTIC® compressive garments are machine or manually washable at the temperature of 30°C max. Before washing, fasten the zip-fasteners and Velcro fasteners. Do not squeeze and iron, dry in a shady and airy place. Drying at a high temperature reduces the elasticity of the material and reduces its curative effect.

DESIGN:



CONTENT

LIST OF THE PRODUCTS

FEMALE VESTS & BODY SUIT

MT smooth	6
MT formed	6
MH special	7
MH	7
MHG	7
MHF	8
MHD	8
MHB	8
MGG	9
MGF	9
MGD	10
MGB	10

FEMALE GIRDLES

VH	11
VH special	11
VH variant	11
VG	12
VG variant	12
VF	13
VF variant	13
VD	14
VD special	14
VD variant	14
VB	15
VB variant	15
TG	16
TG variant	16
TF	17
TF variant	17
TD	18
TD special	18
TD variant	18
TB	19
TB variant	19

BREAST BANDAGES & STABILIZER BANDS

SI	20
SI special	20
SI formed	20
PI formed	21
PS formed	21
PI extra	22
PS extra	22
PI special	23
PS special	23
PI standard	24
PI medical	24
PI super	25
PS variant	25

FACIAL BANDAGES

FM	26
FM formed	26
FM special	26
FM extra	27
FM whole face	27

OTHERS

PU	27
AP	28
AS	28
SK	29
FD	29
KP	30
KP special	30

MALE BANDAGES

MTm	31
MTmL	31
VHm	32
MGm	32



MT smooth



Use: Liposuction of arms, back and chest

Types: with zipper or Comfort

Colours: natural, black



MH special



Use: Mammoplasty, Abdominoplasty, Liposuction of abdomen and back

Types: with zipper or Comfort

Colours: natural, black



MH



Use: Mammoplasty, Abdominoplasty, Liposuction of abdomen, back and upper arms

Types: with zipper or Comfort

Colours: natural, black



MT formed



Use: Liposuction of back, Mammoplasty /Augmentation, Breast reduction and Mastopexy/

Types: with zipper or Comfort

Colours: natural, black



MHG



Use: Mammoplasty, Abdominoplasty, Liposuction of abdomen, back, upper arms, buttocks and hips

Types: with zipper or Comfort

Colours: natural, black



Size	Measurements in cm			
	Breasts circumference	Under breasts	Waist	Arm
SS	80 - 86	68 - 72	58 - 63	26 - 28
S	87 - 92	73 - 77	64 - 69	29 - 30
M	93 - 97	78 - 82	70 - 76	31 - 32
L	98 - 102	83 - 87	77 - 85	33 - 35
XL	103 - 107	88 - 94	86 - 94	36 - 37

Size	Measurements in cm				
	Breasts circumference	Under breasts	Waist	Hips	Thigh
SS	80 - 86	68 - 72	58 - 63	84 - 88	48 - 50
S	87 - 92	73 - 77	64 - 69	89 - 94	51 - 55
M	93 - 97	78 - 82	70 - 76	95 - 101	56 - 60
L	98 - 102	83 - 87	77 - 85	102 - 110	61 - 65
XL	103 - 107	88 - 94	86 - 94	111 - 119	66 - 70



MHF



Use: Mammoplasty, Abdominoplasty, Liposuction of abdomen, back, upper arms, buttocks, hips and thighs

Types: with zipper or Comfort

Colours: natural, black



MGG



Use: Mammoplasty, Abdominoplasty, Liposuction of abdomen, back, buttocks and hips

Types: with zipper or Comfort

Colours: natural, black



MHD



Use: Mammoplasty, Abdominoplasty, Liposuction of abdomen, back, upper arms, buttocks, hips and knees

Types: with zipper or Comfort

Colours: natural, black

Bottom leg part: with zippers or without zippers



MGF



Use: Mammoplasty, Abdominoplasty, Liposuction of abdomen, back, buttocks, hips and thighs

Types: with zipper or Comfort

Colours: natural, black



MHB



Use: Mammoplasty, Abdominoplasty, Liposuction of abdomen, back, upper arms, buttocks, hips, knees and calves

Types: with zipper or Comfort

Colours: natural, black

Bottom leg part: with zippers or without zippers



Size	Measurements in cm						
	Breasts circumference	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	80 - 86	68 - 72	58 - 63	84 - 88	48 - 50	32 - 34	27 - 29
S	87 - 92	73 - 77	64 - 69	89 - 94	51 - 55	35 - 39	30 - 32
M	93 - 97	78 - 82	70 - 76	95 - 101	56 - 60	40 - 44	33 - 35
L	98 - 102	83 - 87	77 - 85	102 - 110	61 - 65	45 - 49	36 - 38
XL	103 - 107	88 - 94	86 - 94	111 - 119	66 - 70	50 - 54	39 - 41

Size	Measurements in cm				
	Breasts circumference	Under breasts	Waist	Hips	Thigh
SS	80 - 86	68 - 72	58 - 63	84 - 88	48 - 50
S	87 - 92	73 - 77	64 - 69	89 - 94	51 - 55
M	93 - 97	78 - 82	70 - 76	95 - 101	56 - 60
L	98 - 102	83 - 87	77 - 85	102 - 110	61 - 65
XL	103 - 107	88 - 94	86 - 94	111 - 119	66 - 70



MGD



Use: Mammoplasty, Abdominoplasty, Liposuction of abdomen, back, buttocks, hips, thighs and knees

Types: with zipper or Comfort

Colours: natural, black

Bottom leg part: with zippers or without zippers



MGB



Use: Mammoplasty, Abdominoplasty, Liposuction of abdomen, back, buttocks, hips, thighs, knees and calves

Types: with zipper or Comfort

Colours: natural, black

Bottom leg part: with zippers or without zippers



Size	Measurements in cm						
	Breasts circumference	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	80 - 86	68 - 72	58 - 63	84 - 88	48 - 50	32 - 34	27 - 29
S	87 - 92	73 - 77	64 - 69	89 - 94	51 - 55	35 - 39	30 - 32
M	93 - 97	78 - 82	70 - 76	95 - 101	56 - 60	40 - 44	33 - 35
L	98 - 102	83 - 87	77 - 85	102 - 110	61 - 65	45 - 49	36 - 38
XL	103 - 107	88 - 94	86 - 94	111 - 119	66 - 70	50 - 54	39 - 41
XXL	108 - 112	93 - 99	91 - 100	120 - 129	71 - 75	55 - 60	42 - 44



VH



Use: Abdominoplasty, Liposuction of abdomen

Types: with zipper, without zipper or Comfort

Colours: natural, black



VH special



Use: Abdominoplasty, Liposuction of abdomen

Types: with zippers, without zippers or Comfort

Colours: natural, black



VH variant



Use: Abdominoplasty, Liposuction of abdomen

Type: variant

Colours: natural, black



Size	SS	S	M	L	XL	XXL
Under breasts	61-66	67-72	73-79	80-88	89-97	98-106
Waist	58-63	64-69	70-76	77-85	86-94	95-103
Hips	81-88	89-94	95-101	102-110	111-119	120-129



VG



Use: Liposuction of abdomen, buttocks, hips and thighs

Types: with zippers, without zippers or Comfort

Colours: natural, black



VG variant



Use: Liposuction of abdomen, buttocks, hips and thighs

Type: variant

Colours: natural, black



Size	Measurements in cm			
	Under breasts	Waist	Hips	Thigh
SS	61 - 66	58 - 63	84 - 88	48 - 50
S	67 - 72	64 - 69	89 - 94	51 - 55
M	73 - 79	70 - 76	95 - 101	56 - 60
L	80 - 88	77 - 85	102 - 110	61 - 65
XL	89 - 97	86 - 94	111 - 119	66 - 70
XXL	98 - 106	95 - 103	120 - 128	71 - 75



VF



Use: Liposuction of abdomen, buttocks, hips and thighs

Types: with zippers, without zippers or Comfort

Colours: natural, black



VF variant



Use: Liposuction of abdomen, buttocks, hips and thighs

Type: variant

Colours: natural, black



Size	Measurements in cm				
	Under breasts	Waist	Hips	Thigh	Above knee
SS	61 - 66	58 - 63	84 - 88	48 - 50	32 - 34
S	67 - 72	64 - 69	89 - 94	51 - 55	35 - 39
M	73 - 79	70 - 76	95 - 101	56 - 60	40 - 44
L	80 - 88	77 - 85	102 - 110	61 - 65	45 - 49
XL	89 - 97	86 - 94	111 - 119	66 - 70	50 - 54
XXL	98 - 106	95 - 103	120 - 128	71 - 75	55 - 59



VD



Use: Liposuction of abdomen, buttocks, hips, thighs and knees

Types: with zippers, without zippers or Comfort

Colours: natural, black



VB



Use: Liposuction of abdomen, buttocks, hips, thighs, knees and calves. Stabilization of calf implants.

Types: with zippers, without zippers or Comfort

Colours: natural, black



VD special



Use: Liposuction of abdomen, buttocks, hips, thighs and knees

Types: with zippers, without zippers or Comfort

Colours: natural, black

Double underlay inside and outside of thighs.



VB variant



Use: Liposuction of abdomen, buttocks, hips, thighs, knees and calves

Type: variant

Colours: natural, black



VD variant



Use: Liposuction of abdomen, buttocks, hips, thighs and knees

Type: variant

Colours: natural, black



Size	Measurements in cm					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	61 - 66	58 - 63	84 - 88	48 - 50	32 - 34	27 - 29
S	67 - 72	64 - 69	89 - 94	51 - 55	35 - 39	30 - 32
M	73 - 79	70 - 76	95 - 101	56 - 60	40 - 44	33 - 35
L	80 - 88	77 - 85	102 - 110	61 - 65	45 - 49	36 - 38
XL	89 - 97	86 - 94	111 - 119	66 - 70	50 - 54	39 - 41
XXL	98 - 106	95 - 103	120 - 128	71 - 75	55 - 60	42 - 44

Size	Measurements in cm					
	Under breasts	Waist	Hips	Thigh	Above knee	Calf
SS	61 - 66	58 - 63	84 - 88	48 - 50	32 - 34	27 - 29
S	67 - 72	64 - 69	89 - 94	51 - 55	35 - 39	30 - 32
M	73 - 79	70 - 76	95 - 101	56 - 60	40 - 44	33 - 35
L	80 - 88	77 - 85	102 - 110	61 - 65	45 - 49	36 - 38
XL	89 - 97	86 - 94	111 - 119	66 - 70	50 - 54	39 - 41
XXL	98 - 106	95 - 103	120 - 128	71 - 75	55 - 60	42 - 44



TG



Use: Liposuction of abdomen, buttocks, hips, thighs

Types: with zippers, without zippers or Comfort

Colours: natural, black



TF



Use: Liposuction of abdomen, buttocks, hips, thighs

Types: with zippers, without zippers or Comfort

Colours: natural, black



TG variant



Use: Liposuction of abdomen, buttocks, hips, thighs

Type: variant

Colours: natural, black



TF variant



Use: Liposuction of abdomen, buttocks, hips, thighs

Type: variant

Colours: natural, black



Size	Measurements in cm		
	Waist	Hips	Thigh
SS	58 - 63	84 - 88	48 - 50
S	64 - 69	89 - 94	51 - 55
M	70 - 76	95 - 101	56 - 60
L	77 - 85	102 - 110	61 - 65
XL	86 - 94	111 - 119	66 - 70
XXL	95 - 103	120 - 128	71 - 75

Size	Measurements in cm			
	Waist	Hips	Thigh	Above knee
SS	58 - 63	84 - 88	48 - 50	32 - 34
S	64 - 69	89 - 94	51 - 55	35 - 39
M	70 - 76	95 - 101	56 - 60	40 - 44
L	77 - 85	102 - 110	61 - 65	45 - 49
XL	86 - 94	111 - 119	66 - 70	50 - 54
XXL	95 - 103	120 - 128	71 - 75	55 - 59



TD



Use: Liposuction of abdomen, buttocks, hips, thighs and knees

Types: with zippers, without zippers or Comfort

Colours: natural, black



TB



Use: Liposuction of abdomen, buttocks, hips, thighs, knees and calves. Stabilization of calf implants.

Types: with zippers, without zippers or Comfort

Colours: natural, black



TD special



Use: Liposuction of abdomen, buttocks, hips, thighs and knees

Types: with zippers, without zippers or Comfort

Colours: natural, black

Double underlay inside and outside of thighs.



TB variant



Use: Liposuction of abdomen, buttocks, hips, thighs, knees and calves. Stabilization of calf implants.

Type: variant

Colours: natural, black



TD variant



Use: Liposuction of abdomen, buttocks, hips, thighs and knees

Type: variant

Colours: natural, black



Size	Measurements in cm				
	Waist	Hips	Thigh	Above knee	Calf
SS	58 - 63	84 - 88	48 - 50	32 - 34	27 - 29
S	64 - 69	89 - 94	51 - 55	35 - 39	30 - 32
M	70 - 76	95 - 101	56 - 60	40 - 44	33 - 35
L	77 - 85	102 - 110	61 - 65	45 - 49	36 - 38
XL	86 - 94	111 - 119	66 - 70	50 - 54	39 - 41
XXL	95 - 103	121 - 129	71 - 75	55 - 59	42 - 44

Size	Measurements in cm				
	Waist	Hips	Thigh	Above knee	Calf
SS	58 - 63	84 - 88	48 - 50	32 - 34	27 - 29
S	64 - 69	89 - 94	51 - 55	35 - 39	30 - 32
M	70 - 76	95 - 101	56 - 60	40 - 44	33 - 35
L	77 - 85	102 - 110	61 - 65	45 - 49	36 - 38
XL	86 - 94	111 - 119	66 - 70	50 - 54	39 - 41
XXL	95 - 103	121 - 129	71 - 75	55 - 59	42 - 44



SI

Use: Implant Stabilizer Band for a separate use or in combination with surgical bras type PI as a stabilization support of implants following Augmentation

Colours: white, black

Size	Breasts circumference in cm
S	68 - 77
M	78 - 87
L	88 - 97
XL	98 - 105



SI special

Use: Implant Stabilizer Band for a separate use or in combination with surgical bras type PI as a stabilization support of implants following Augmentation

Colour: white

Size	Breasts circumference in cm
S	68 - 77
M	78 - 87
L	88 - 97
XL	98 - 105



SI formed

Use: Implant Stabilizer Band for a separate use or in combination with surgical bras type PI as a stabilization support of implants following Augmentation

Colours: white, black

Size	Breasts circumference in cm	Under breasts in cm
60	72 - 80	58 - 62
65	77 - 85	63 - 67
70	82 - 90	68 - 72
75	87 - 95	73 - 77
80	92 - 100	78 - 82
85	97 - 105	83 - 87
90	102 - 110	88 - 92
95	107 - 115	93 - 97
100	112 - 120	98 - 102



PI formed



Use: Mammoplasty /Augmentation, Breast reduction or Mastopexy/
Colours: white, black



PS formed



Use: Augmentation
Colours: white, black



Size	Under breasts	Breasts circumference in cm			
		Cup A	Cup B	Cup C	Cup D
60	58 - 62	72 - 74	74 - 76	76 - 78	78 - 80
65	63 - 67	77 - 79	79 - 81	81 - 83	83 - 85
70	68 - 72	82 - 84	84 - 86	86 - 88	88 - 90
75	73 - 77	87 - 89	89 - 91	91 - 93	93 - 95
80	78 - 82	92 - 94	94 - 96	96 - 98	98 - 100
85	83 - 87	97 - 99	99 - 101	101 - 103	103 - 105
90	88 - 92	102 - 104	104 - 106	106 - 108	108 - 110
95	93 - 97	107 - 109	109 - 111	111 - 113	113 - 115
100	98 - 102	112 - 114	114 - 116	116 - 118	118 - 120



PI extra

Use: Mammoplasty /Augmentation,
Breast reduction or Mastopexy/
Colour: white



PS extra

Use: Augmentation
Colour: white



PI special

Use: Mammoplasty /Augmentation,
Breast reduction or Mastopexy/
Colours: white, black



PS special

Use: Augmentation
Colours: white, black



Size	Measurements in cm	
	Breasts circumference	Under breasts
SS	80 - 86	68 - 72
S	87 - 92	73 - 77
M	93 - 97	78 - 82
L	98 - 102	83 - 87
XL	103 - 107	88 - 92
XXL	108 - 112	93 - 97

Size	Under breasts	Breasts circumference in cm			
		Cup A	Cup B	Cup C	Cup D
60	58 - 62	72 - 74	74 - 76	76 - 78	78 - 80
65	63 - 67	77 - 79	79 - 81	81 - 83	83 - 85
70	68 - 72	82 - 84	84 - 86	86 - 88	88 - 90
75	73 - 77	87 - 89	89 - 91	91 - 93	93 - 95
80	78 - 82	92 - 94	94 - 96	96 - 98	98 - 100
85	83 - 87	97 - 99	99 - 101	101 - 103	103 - 105
90	88 - 92	102 - 104	104 - 106	106 - 108	108 - 110
95	93 - 97	107 - 109	109 - 111	111 - 113	113 - 115
100	98 - 102	112 - 114	114 - 116	116 - 118	118 - 120



PI standard



Use: Mammoplasty /Augmentation,
Breast reduction or Mastopexy/
Colours: white, black, natural



Size	Under breasts	Breasts circumference in cm				
		Cup AA	Cup A	Cup B	Cup C	Cup D
60	58 - 62	70 - 72	72 - 74	74 - 76	76 - 78	78 - 80
65	63 - 67	75 - 77	77 - 79	79 - 81	81 - 83	83 - 85
70	68 - 72	80 - 82	82 - 84	84 - 86	86 - 88	88 - 90
75	73 - 77	85 - 87	87 - 89	89 - 91	91 - 93	93 - 95
80	78 - 82	-	92 - 94	94 - 96	96 - 98	98 - 100
85	83 - 87	-	97 - 99	99 - 101	101 - 103	103 - 105
90	88 - 92	-	102 - 104	104 - 106	106 - 108	108 - 110
95	93 - 97	-	107 - 109	109 - 111	111 - 113	113 - 115
100	98 - 102	-	112 - 114	114 - 116	116 - 118	118 - 120



PI super

Use: Mammoplasty /Augmentation,
Breast reduction or Mastopexy/
Colours: white, black



Size	Under breasts	Breasts circumference in cm			
		Cup A	Cup B	Cup C	Cup D
60	58 - 62	72 - 74	74 - 76	76 - 78	78 - 80
65	63 - 67	77 - 79	79 - 81	81 - 83	83 - 85
70	68 - 72	82 - 84	84 - 86	86 - 88	88 - 90
75	73 - 77	87 - 89	89 - 91	91 - 93	93 - 95
80	78 - 82	92 - 94	94 - 96	96 - 98	98 - 100
85	83 - 87	97 - 99	99 - 101	101 - 103	103 - 105
90	88 - 92	102 - 104	104 - 106	106 - 108	108 - 110
95	93 - 97	107 - 109	109 - 111	111 - 113	113 - 115
100	98 - 102	112 - 114	114 - 116	116 - 118	118 - 120



PI medical

Use: Mammoplasty /Augmentation,
Breast reduction or Mastopexy/
Colours: white, black



Size	Under breasts	Breasts circumference in cm				
		Cup AA	Cup A	Cup B	Cup C	Cup D
60	58 - 62	70 - 72	72 - 74	74 - 76	76 - 78	78 - 80
65	63 - 67	75 - 77	77 - 79	79 - 81	81 - 83	83 - 85
70	68 - 72	80 - 82	82 - 84	84 - 86	86 - 88	88 - 90
75	73 - 77	85 - 87	87 - 89	89 - 91	91 - 93	93 - 95
80	78 - 82	-	92 - 94	94 - 96	96 - 98	98 - 100
85	83 - 87	-	97 - 99	99 - 101	101 - 103	103 - 105
90	88 - 92	-	102 - 104	104 - 106	106 - 108	108 - 110
95	93 - 97	-	107 - 109	109 - 111	111 - 113	113 - 115
100	98 - 102	-	112 - 114	114 - 116	116 - 118	118 - 120



PS variant

Use: Augmentation
Colour: white



Size	Measurements in cm	
	Breasts circumference	Under breasts
SS	80 - 86	68 - 72
S	87 - 92	73 - 77
M	93 - 97	78 - 82
L	98 - 102	83 - 87
XL	103 - 107	88 - 92
XXL	112 - 116	93 - 97



FM



Use: Face, neck and chin operations,
chin implants surgery

Colours: natural, black



FM extra

Use: Face, neck and chin operations,
chin implants surgery

Colour: white

Size	Head measurement in cm /chin - top of the head/
SS	50 - 53
S	54 - 59
M	60 - 65
L	66 - 70
XL	71 - 74



FM formed



Use: Face, neck and chin operations,
chin implants surgery

Colours: natural, black



FM whole face

Use: Face, neck and chin operations,
chin implants surgery, burns

Colours: natural, black

Size	Head measurement in cm /chin - top of the head/
SS	50 - 53
S	54 - 59
M	60 - 65
L	66 - 70
XL	71 - 74



FM special



Use: Face, neck and chin operations,
chin implants surgery, ear reconstruction

Colours: natural, black



PU

Use: Otoplasty, ear reconstruction

Colours: natural, black

Size	Head measurement in cm /over ears/
SS	50 - 52
S	53 - 55
M	56 - 58
L	59 - 61
XL	63 - 64



Size	Head measurement in cm /chin - top of the head/
SS	50 - 53
S	54 - 59
M	60 - 65
L	66 - 70
XL	71 - 74



AP

Use: Liposuction of armpit and arms
Colours: natural, black



SK

Use: Mould the figure, making it firmer and adding the tone, hold shape perfectly, make the tummy flat and shape up the hips
Colours: natural, black, white



Size	Confection size	Waist in cm	Hips in cm
SS	32 - 34	58 - 63	84 - 88
S	36 - 38	64 - 69	89 - 94
M	40 - 42	70 - 76	95 - 101
L	44 - 46	77 - 85	102 - 110
XL	48 - 50	86 - 94	111 - 119
XXL	52 - 54	95 - 103	120 - 128



AS

Use: Liposuction of arms
Colours: natural, black



Size	Measurements in cm	
	Arm circumference	Under elbow circumference
SS	26 - 28	21 - 23
S	29 - 30	24 - 25
M	31 - 32	26 - 27
L	33 - 35	28 - 29
XL	36 - 37	30 - 31
XXL	39 - 39	33 - 33



FD

Use: Liposuction of knees
Colours: natural, black



Size	Measurements in cm	
	Above knee	Below knee
SS	32 - 34	27 - 29
S	35 - 39	30 - 32
M	40 - 44	33 - 35
L	45 - 49	36 - 38
XL	50 - 54	39 - 41
XXL	55 - 59	42 - 44



KP

Use: Abdominoplasty
Colours: white, black



MTm

Use: Plastic surgeries of back and chest
Types: with zipper or Comfort
Colours: natural, black



KP special

Use: Abdominoplasty
Colour: white
Type: side reinforcement by flexible strap



MTmL

Use: Plastic surgeries of back, chest and shoulders
Types: with zipper or Comfort
Colours: natural, black



Size	Waist in cm	Binder height in cm
SS	58 - 63	16, 23, 31, 39
S	64 - 69	16, 23, 31, 39
M	70 - 76	16, 23, 31, 39
L	77 - 85	16, 23, 31, 39
XL	86 - 94	16, 23, 31, 39

Size	Measurements in cm	
	Chest	Waist
SS	87 - 89	75 - 77
S	90 - 93	78 - 81
M	94 - 101	82 - 89
L	102 - 109	90 - 99
XL	110 - 117	100 - 109



VHm



Use: Abdominoplasty, Liposuction of abdomen and hips

Types: with zipper, without zipper or Comfort

Colours: natural, black

Size	Measurements in cm					
	SS	S	M	L	XL	XXL
Under breasts	87-89	90-93	94-101	102-109	110-117	118-125
Waist	75-77	78-81	82-89	90-99	100-109	110-119
Hips	82-85	86-89	90-96	97-104	105-112	113-120



MGm

Use: Liposuction of abdomen, chest, back, shoulders and hips

Types: with zipper or Comfort

Colours: natural, black

Size	Measurements in cm			
	Chest	Waist	Hips	Thigh
SS	87 - 89	75 - 77	82 - 85	47 - 50
S	90 - 93	78 - 81	86 - 89	51 - 55
M	94 - 101	82 - 89	90 - 96	56 - 60
L	102 - 109	90 - 99	97 - 104	61 - 65
XL	110 - 117	100 - 109	105 - 112	66 - 68
XXL	118 - 125	110 - 119	113 - 120	69 - 71



lipoelastic®

LIPOELASTIC a. s.
Krhová 53
757 01 Valašské Meziříčí
CZECH REPUBLIC
tel: +420 571 633 512
fax: +420 571 633 268
e-mail: info@lipoelastic.com
<http://www.lipoelastic.com>



October 2008